

Instructions: MRI Rating Tool

Project Description

The simple MRI rating tool allows users to open MRI images and provide a numerical rating that describes the impact of artifact on the images quality.

The instructions below detail how to use this rating tool to assess MRI NIfTI files.

Requirements:

- ITK-SNAP (medical image viewing software)
 - <http://www.itksnap.org/pmwiki/pmwiki.php?n=Downloads.SNAP3>
- Python (latest stable version recommended)
 - <https://www.python.org/downloads/>

SECTION 1: FILE ORGANISATION

Ensure all necessary files are in the correct directory


For the MRI rating tool to run correctly the following items must be present in the same directory

- Python file containing the objects for the simpledb database management tool (simpledb_objects.py)
- Python file containing the functional rating tool script (mri_qa-marker.py)
- Directory of images in NIfTI file format

Step 1: Download python files (.py)

- Create a new directory (i.e. folder) to store files for the rating tool
- Download the following files and place them into the directory:
 - mri_qa-marker.py
 - simpledb_objects.py

Image: In the example below they have been moved into a directory named 'code'




Step 2: Move NIfTI files into this directory

- Move the images (NIfTI files) into the directory containing the two python scripts

Note: the directory containing the images can have any name, the script will find any NIfTI (nii.gz) files contained in directory where the python scripts are stored

Image: In the example below a directory named 'my_mri_images' was moved into the directory



SECTION 2: RUN THE SCRIPT

To run the script you will need to have python installed, if you are unsure whether you have python installed simply follow step 1 below. If it is installed the script will run correctly, otherwise it will return the error: **'python: command not found'** or **'python3: command not found'**. For details about how to install python please see: <https://www.python.org/downloads/>

Step 1: Open a Terminal in the directory containing the python files and the MRI images

- Open a terminal and ensure you are in the directory containing the files
 - **Example:** <your_username>:~/your/file/path/code\$
 - In this example the terminal is open in the 'code' directory where the files are contained
- Execute the script by typing into the terminal and pressing enter:
 - `python mri_qa-marker.py`

Note: if you are using python version 3 you will need to type: **python3 mri_qa-marker.py**

Step 2: Choose your session

- You will be prompted to enter whether you would like to resume a previous review session or start a new one, indicate your choice by entering the corresponding number

Note: if you have not previously used the rating tool, a directory titled 'mri_rating_record' (see image below) will be generated, this is the output directory where your ratings will be stored as a .json file

```
mollyi@95201-004726-1:~/Documents/Projects/SimpleQA/code$ python3 mri_qa-marker.py


Directory 'MRI_rating_record' created: your ratings will be stored here.

Would you like to resume a previous session?
1 - New session
2 - Resume previous sessions
Selection: 1
```

SECTION 3A: CHOOSING A NEW SESSION

Select a new session by entering '1'. This will bring up the ITK-SNAP window and your rating should be entered into the terminal.

Image: ITK-SNAP will open in a separate window (left), your rating should be entered in the terminal (right)




Step 1: Provide your rating

- Input your rating, a numerical rating between 1 and 5, into the terminal and press the 'Enter' key
- After pressing 'Enter' your rating will displayed back to you in the terminal (See image below)
- The current image will be closed and the next image for review will automatically be opened
- Continue through the images providing your rating in the terminal

Note: there is no need to close the ITK-SNAP window, it will close automatically when you provide your rating and the next image to be reviewed will be opened

Image: your rating will be displayed back to you in the terminal ('Your rating was 4')



SECTION 3B: RESUMING A PREVIOUS SESSION

Resume reviewing a previous session by entering '2'. You will be prompted to select the session you wish to resume. The ITK-SNAP window will then only open images that have not yet been reviewed.

Step 1: Select your session

- Choose from the list of existing rating files by entering the corresponding number into the terminal


Image: There are 2 available sessions to resume ('MRIrate_session-25-03-2022_16:23:03' and 'MRIrate_session-25-03-2022_16:13:54')

```
mollyi@9520l-004726-l: ~/Documents/Projects/SimpleQA/code
mollyi@9520l-004726-l:~/Documents/Projects/SimpleQA/code$ python3 mri_qa-marker.py
Would you like to resume a previous session?
1 - New session
2 - Resume previous sessions
Selection: 2
1. MRIrate_session-25-03-2022_16:23:03
2. MRIrate_session-25-03-2022_16:13:54
Enter the number of the session you wish to resume: █
```

Step 2: Provide your rating

- After selecting the session to resume, the images that have been previously reviewed will appear in the terminal with the comment "file <file_name> has already been rated" (see image below)
- Once the script encounters an image file that has not yet been rated, you will be prompted to enter your rating and ITK-SNAP window containing the file will open
- Input your rating, a numerical rating between 1 and 5, into the terminal and press the 'Enter' key
- The previous image will be closed and the next image for review will automatically be opened
- Continue through the images providing your rating in the terminal

Note: there is no need to close the ITK-SNAP window, it will close automatically when you provide your rating and the next image to be reviewed will be opened



SECTION 4: ENDING YOUR SESSION


- There are 2 ways to end the session:
 1. Review all images in the directory
 2. Typing 'Ctrl' + 'C' in the terminal during a review session

Note: the files you review will be stored even if you terminate your session early

- In both cases the text '--END OF IMAGES--' will appear and you will be given an option to convert your review file to a .csv format (see image below)
 - If you select 'y' to convert the files, your .csv file can be found in the 'MRI_rating_record' directory with the accompanying .json file

Note: it is not necessary to make a .csv file to save your review, it is automatically saved as a .json file

Image: See example below where 'Ctrl' + 'C' was entered into the terminal to end the session

A terminal window with a dark purple background and light green text. The text shows a prompt 'Rating (1 to 5) : ^C', followed by '--END OF IMAGES--', then 'Output the .JSON file to a .CSV file?', and finally 'Press \'y\' to convert or Enter to quit: ' with a cursor at the end.

```
Rating (1 to 5) : ^C
--END OF IMAGES--
Output the .JSON file to a .CSV file?
Press 'y' to convert or Enter to quit: 
```

SECTION 5: ACCESSING YOUR REVIEW FILES

- All review files (.json, .csv) can be found in the directory 'MRI_rating_record' contained within the directory where the python scripts can be found